



Assembly Directions

- 1) Cut on all of the solid line.
- 2) Make certain to cut the solid line in section C (this should be a small slit in the paper which light can enter spectroscope) and the box in section F
- 3) Attach the piece of CD to the spot labelled.
- 4) Fold on all of the dotted lines
- 5) Fold flaps A and B up at a 90 degree angle.
- 6) Fold sections C, D, E, and F along the edges of flaps A and B.
- 7) Fold the small flaps attached to each section over to cover any holes from which light could enter.
- 8) Look through the hole in section F at the CD fragment. You should be able to see a rainbow from the light being split by the grating on the CD.

